

# alma MAMMO

## MAMMOGRAPHY SCREENING

Alma MAMMO is specifically designed for the display and evaluation of mammography screening. This module automatically classifies images according to a projection (cranial-caudal, half lateral oblique, lateral, etc.), and laterally (left or right) and displays them into an organized manner on the monitors.

Alma MAMMO allows each user to create their own workflow, to suit the image review protocol established by each user or centre.

Alma MAMMO also allows creating your own workflow to suit the best image review protocol established by each user. These tools will allow increased speed and accuracy for the diagnosis and review.



The product offers a powerful DICOM study manager that allows not only to query and retrieve studies from a PACS server, but also to send new series generated during the evaluation of the study.

The product is implemented according to recommendations of IHE international guidelines, and obtained accreditation as Image Display at the Connectathon organized by IHE Europe.

**Configuration of workflow.**

**Exploration of previous images.**

**CAD mark readings.**

**Interactive magnifying glass.**

**Hanging protocol configuration.**