


MULTIPLANAR RECONSTRUCTION

Alma MPR is a multiplanar reconstruction module. Based on a sequence of images, it constructs a three-dimensional model and allows the simultaneous visualisation of MPR projections in Axial, Coronal, Sagittal and Oblique plane in a completely interactive way.

The visualisation of oblique MPR projections allows us to obtain three-dimensional sections of the model in whichever plane you prefer.



The software allows the client to interactively modify the plane of the section and obtain the corresponding image in real time.

Alma MPR, validated by the IHE, is totally compatible with the DICOM standard and is completely prepared to work integrally with HIS/RIS/PACS clinical information systems.

- ★ Zoom, Pann and W/L.
- ★ Creation of new series.
- ★ Annotations.
- ★ Oblique MPR.
- ★ Slice thickness (MIP, MinIP, VR, Avg).
- ★ Curved MPR.
- ★ Library of palettes.